


All in BeamCheck

- More efficient workflows
- More precise analysis results
- More flexible measurement settings




BeamCheck supports the following procedures:

Procedures

Center	Geometry	Dose Characteristics	Evaluation
<ul style="list-style-type: none"> ✓ Point Dose ✓ Field Center Deviation X/Y ✓ Wedge Factor 	<ul style="list-style-type: none"> ✓ Size ✓ Field Edge -/+ ✓ Penumbra -/+ ✓ Central Leaf Position ✓ Leaf Average Position -/+ 	<ul style="list-style-type: none"> ✓ Flatness ✓ Symmetry ✓ PDD₁₀ ✓ PDD₂₀ ✓ R₅₀ 	<ul style="list-style-type: none"> ✓ Mean Deviation ✓ Gamma ✓ Leaf Position Pass Rate




BeamCheck Wedge



BeamCheck Fixture

BeamCheck is suitable for monthly checks of LINAC and also helps users with mechanical inspections and post-repair commissioning.

Ion Chamber Dual-Layer Arrangement



Ion Chamber Quantity: 759
Field Size: 30cm x 30cm
Max Axis Range (x/y): 31.4cm x 31.4cm


X-axis(Dual-layer arrangement): 157 ion chambers;
Y-axis(Dual-layer arrangement): 150 ion chambers.
Used for measuring output, flatness, symmetry, etc.
Ion chamber spacing = 2mm(center to center)

Four columns(Y1,Y2,Y3,Y4;Dual-layer arrangement):
Each column has 104 ion chambers:
Used for measuring leaf positioning accuracy.
Ion chamber spacing=2.5mm
(center to center)
MLC measurement accuracy=1mm

10cm x 10cm & 20cm x 20cm
Total of 36 ion chambers on the boundary:
Used for checking coincidence of light field and radiation field.

Built-in battery, no external power supply required during measurement.

Charges via Type-C interface, providing up to 9 hours of continuous usage when fully charged.



XBeam: BeamCheck software

Toolbar
One-click operation, convenient and efficient.

Measurement Conditions
Set anytime, start measurement freely.

Ion Chamber Measurement Curve
Real-time monitoring, intuitive display.

Analysis Result Area

Item	Value	Unit
尺寸 [mm]	10.04	-0.20
辐射剂量率 \dot{D} [mSv/h]	-5.03/5.01	0.10%/0.10
半影 r_{50} [cm]	0.55/0.56	-0.10/0.00
中心轴射线位置 r_{50} [cm]	-	-
叶片半轴位置 r_{50} [cm]	-	-
平坦度 [%]	1.29	0.03
对称性 [%]	-14.43	0.01
PDD [%]	-	-
R50 [%]	-	-
平均偏置 [%]	-	0.03

Four Core Functions

- Project-Based Management**
 Measurement tasks are categorized by time, ensuring clearer data management and easier tracking.
- Real-Time Analysis**
 Instantly generate measurement curves for an intuitive display, aiding quick decision-making.
- Custom Measurement Settings**
 Flexibly adjust measurement conditions to suit different equipment and QA processes.
- Comprehensive Reports**
 Provides detailed measurement data analysis and reports covering 17 key indicators.

ID	Name	Machine	Volume	Spacing	Date	Status
1	#####	3080000	GMV	10x10	23-05-16 13:40	运行
2	#####	3080000	GMV	10x10	23-05-16 13:40	运行
3	#####	3080000	GMV	10x10	23-05-16 13:40	运行
4	#####	3080000	GMV	10x10	23-05-16 13:40	运行
5	#####	3080000	GMV	10x10	23-05-16 13:40	运行
6	#####	3080000	GMV	10x10	23-05-16 13:40	运行
7	#####	3080000	GMV	10x10	23-05-16 13:40	运行
8	#####	3080000	GMV	10x10	23-05-16 13:40	运行
9	#####	3080000	GMV	10x10	23-05-16 13:40	运行

Technical Parameters

Detector Type Ion chamber	Detector Quantity 759	Detector Volume 0.051cm³	Detector Spacing 2mm(Min)
Nominal Sensitivity 2.2nC/Gy±0.2nC/Gy	Field Size 30cm×30cm	Repeatability ±0.5%	Non-linearity ±0.5%
Dimensions (L×W×H) 435mm×341mm×22mm	Weight 5kg		